



Premium micropipette tips from Xcelris

- **Uncompromised Production:** Manufactured in fully automatic, touch free clean room production facility
- **High Quality Plastic:** Produced from medical grade virgin polypropylene plastic
- **Contamination free:** All tips are free from DNase, RNase, Pyrogen and PCR inhibitors
- **Compatible Design:** Innovative design compatible with pipettes from other manufacturers
- **Wide variety of Packaging Options:** Available in Bulk, Rack and Reloading Stack
- **High Transparency with Readable marks:** Ensure error free pipetting with visible graduated marks on micropipette tips

Low Retention Tips

Ensure Maximum Accuracy with Low Retention Tips

Manufactured from specifically developed resins combined with high quality diamond finish mould creates **SUPER HYDROPHOBIC** internal surface



Ordering Information

BULK-Double sealed polyethylene bag with easy tearing

10µl 1000 Tips/Bag (40 Bags/Case)	200µl 1000 Tips/Bag (20 Bags/Case)	1000µl 1000 Tips/Bag (10 Bags/Case)	Low Retention (L)	Colored Tips (C)
XPet-10-B	XPet-200-B	XPet-1000-B	–	–
–	XPet-200-B-C	XPet-1000-B-C	–	✓
XPet-10-L-B	XPet-200-L-B	XPet-1000-L-B	✓	–
–	Xpet-200-L-B-C	Xpet-1000-L-B-C	✓	✓



RACK-Colour coded, fully autoclavable, Designed for single hand operations

10µl 96 Tips x 10 Racks/Pack (10 Packs/Case)	200µl 96 Tips x 10 Racks/Pack (10 Packs/Case)	1000µl 96 Tips x 6 Racks/Pack (10 Packs/Case)	Sterile (S)	Low Retention (L)
XPet-10-S-R	XPet-200-S-R	XPet-1000-S-R	✓	–
XPet-10-L-R	XPet-200-L-R	XPet-1000-L-R	–	✓



RELOADING STACK-Effortless preloading of Racks in prearranged 96-well format

10µl 96 Tips x 10 Plates/Pack (10 Packs/Case)	200µl 96 Tips x 10 Plates/Pack (10 Packs/Case)	1000µl 96 Tips x 5 Plates/Pack (10 Packs/Case)	Sterile (S)	Low Retention (L)
XPet-10-S-K	XPet-200-S-K	XPet-1000-S-K	✓	–
XPet-10-L-K	XPet-200-L-K	XPet-1000-L-K	–	✓

